

Human ARG2 (His tag) recombinant protein



Catalog Number: 519851

General Information

Gene Name Synonym

Arginase-2, mitochondrial; Kidney-type arginase; Non-hepatic arginase; Type II arginase; ARG2

Protein Construction

A DNA sequence encoding the Human ARG2 His24-Gly330 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 10mM TrisHCl, 10mM NaCl, 1mM β -mercaptoethanol, pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ARG2 (His tag) recombinant protein