

LAS1L antibody

Catalog Number: 112151

Product name

LAS1L antibody

IP: 1:500-1:5000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LAS1L Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HEK-293 cells, COLO 320 cells, HeLa cells, HepG2 cells, PC-3 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 83 kDa, 70-75 kDa

Preparation

This antibody was obtained by immunization of LAS1L recombinant protein (Accession Number: NM_031206). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112151(LAS1L antibody) at dilution of 1:1000



IP Result of anti-LAS1L (IP:Catalog No:112151, 4ug; Detection:Catalog No:112151 1:1000) with HEK-293 cells lysate 1200ug.



Immunofluorescent analysis of HepG2 cells using Catalog No:112151(LAS1L Antibody) at dilution of

LAS1L antibody



Catalog Number: 112151

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)