

PPP1R14B antibody



Catalog Number: 114151

Product name

PPP1R14B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPP1R14B Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human heart tissue, mouse liver tissue. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of PPP1R14B recombinant protein (Accession Number: NM_138689). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114151(PPP1R14B antibody) at dilution of 1:400