

Anti-Smad3 Rabbit antibody



Catalog Number: 167551

Product name

Anti-Smad3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Smad3 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Smad3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions


WB: 1/1000

IP: 1/20

Validations



Western blot detection of Smad3 in K562,CHO-K1 cell lysates using Smad3 Rabbit pAb(1:1000 diluted).Predicted band size:48KDa.Observed



band size:52KDa.

Western blot detection of Smad3 in Jurkat cell lysates using Smad3 Rabbit pAb(1:1000 diluted).Predicted band size:48KDa.Observed band size:52KDa.