

# METTL2A Polyclonal Antibody



Catalog Number: 161951

---

## Product name

METTL2A Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to METTL2A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 139-378 of human METTL2A (NP\_859076.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:100

## Validations



Immunofluorescence - METTL2A Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using METTL2A antibody . Blue: DAPI for nuclear staining.