

NARFL antibody



Catalog Number: 112951

Product name

NARFL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NARFL Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 55kDa, 42kDa, 29kDa

Preparation

This antibody was obtained by immunization of NARFL recombinant protein (Accession Number: NM_022493). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112951(NARFL Antibody) at dilution of 1:300