

PLD1 antibody

Catalog Number: 113951

Product name

PLD1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLD1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of PLD1 recombinant protein (Accession Number: XM_005247533). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113951(PLD1 antibody) at dilution of 1:500



IP Result of anti-PLD1 (IP:Catalog No:113951, 4ug; Detection:Catalog No:113951 1:500) with HeLa cells lysate 1200ug.



Immunofluorescent analysis of HeLa cells using Catalog No:113951(PLD1 Antibody) at dilution of

PLD1 antibody



Catalog Number: 113951

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)