

# Anti-ATM Rabbit antibody



Catalog Number: 166951

## Product name

Anti-ATM Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

ATM Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human ATM

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP, FC

## Dilutions


WB: 1/500-1/2000

ICC/IF: 1/50-1/200

IP: 1/80

FC: 1/100

## Validations



Western blot analysis of ATM expression in 293 cell lysate.