

ZNF154 antibody



Catalog Number: 116951

Product name

ZNF154 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF154 Rabbit Polyclonal antibody. Positive WB detected in COS-7 cells, A549 cells, HEK-293 cells, MCF-7 cells. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of ZNF154 recombinant protein (Accession Number: NM_001085384). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COS-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116951(ZNF154 Antibody) at dilution of 1:600