

PSMG2 antibody

Catalog Number: 114251

Product name

PSMG2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PSMG2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer.

Preparation

This antibody was obtained by immunization of PSMG2 recombinant protein (Accession Number: NM_020232). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:114251(PSMG2 antibody) at dilution of 1:50 (under 10x lens)