

RNASE11 Polyclonal Antibody



Catalog Number: 160651

Product name

RNASE11 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to RNASE11

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human RNASE11 (NP_660293.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - RNASE11 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using RNASE11 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - RNASE11 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using RNASE11 antibody at dilution of 1:100 (40x lens).