

PDCD7 antibody



Catalog Number: 113651

Product name

PDCD7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDCD7 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of PDCD7 recombinant protein (Accession Number: BC016992). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113651(PDCD7 Antibody) at dilution of 1:600