

# RGSL1 antibody



Catalog Number: 114651

## Product name

RGSL1 antibody

## Immunogen

[Human RGSL1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RGSL1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 125 kDa

## Preparation

This antibody was obtained by immunization of RGSL1 recombinant protein (Accession Number: NM\_001137669). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114651(RGSL1 Antibody) at dilution of 1:600