

# PPP2R4 Polyclonal Antibody



Catalog Number: 166651

## Product name

PPP2R4 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to PPP2R4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human PPP2R4 (NP\_821067.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations



Immunohistochemistry - PPP2R4 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using PPP2R4 Antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - PPP2R4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney cancer using PPP2R4 Antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - PPP2R4 Polyclonal Antibody


# PPP2R4 Polyclonal Antibody



Catalog Number: 166651

---

Immunohistochemistry of paraffin-embedded human esophageal cancer using PPP2R4 Antibody



at dilution of 1:200 (40x lens).

Immunohistochemistry - PPP2R4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using PPP2R4 Antibody at dilution of 1:200 (40x lens).