

Anti-MROH7 antibody



Catalog Number: 176651

Product name

Anti-MROH7 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to MROH7

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:1,000

IHC-P: 1:50-1:200

Validations




Fig1: Western blot analysis of C1orf175 on recombinant protein using anti-C1orf175 antibody at 1/1,000 dilution.




Fig2: Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-C1orf175 antibody. Counter stained with hematoxylin.