

# EPHA6-Specific antibody



Catalog Number: 110351

## Product name

EPHA6-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EPHA6-Specific Rabbit Polyclonal antibody.  
Positive IHC detected in mouse brain tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001080448).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:110351 (EPHA6-Specific Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:110351 (EPHA6-Specific Antibody) at dilution of 1:50 (under 40x lens)