

Anti-TLE4 antibody

Catalog Number: 176351

Product name

Anti-TLE4 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TLE4

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tle4. aa range:161-210

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using TLE4 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).