

Anti-TLE4 antibody



Catalog Number: 176351

Product name

Rabbit IgG

Anti-TLE4 antibody

Applications

Specificity

WB, IHC-P

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1:500-1:2000


Rabbit polyclonal antibody to TLE4

IHC-P: 1:100-1:300

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tle4. aa range:161-210

Validations



Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Fig1:; Western Blot analysis of various cells using TLE4 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).