

# TRAV20 antibody



Catalog Number: 116351

## Product name

TRAV20 antibody

## Specificity

human; other species not tested.

## Antibody description

TRAV20 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IHC detected in human lymphoma. Observed molecular weight by Western-blot: 31kd

## Preparation

This antibody was obtained by immunization of TRAV20 recombinant protein (Accession Number: BC025727). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116351(TRAV20 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116351(TRAV20 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116351(TRAV20 antibody) at dilution of 1:50 (under 40x lens)