

# Integrin alpha-3 antibody



Catalog Number: 107351

## Product name

Integrin alpha-3 antibody

## Specificity

Human; other species not tested.

## Antibody description

Integrin alpha-3 Mouse Monoclonal antibody. Positive WB detected in PC-3 cells. Positive FC detected in A549 cells. Observed molecular weight by Western-blot: 120 kDa

## Preparation

This antibody was obtained by immunization of Integrin alpha-3 recombinant protein (Accession Number: NM\_002204). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107351(ITGA3 antibody) at dilution of 1:1000



1X10<sup>6</sup> A549 cells were stained with 0.2ug Integrin alpha-3 antibody (Catalog No:107351, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1500.