

Human UBE2V1/CROC1/UBE2V/UEV1/P/OKcl.19 (His tag) recombinant protein



Catalog Number: 518351

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 Variant 1, UEV-1, CROC-1, TRAF6-Regulated IKK Activator 1 Beta Uev1A, UBE2V1, CROC1, UBE2V, UEV1, P/OKcl.19

Protein Construction

A DNA sequence encoding the Human UBE2V1/CROC1/UBE2V/UEV1/P/OKcl.19 Ala2-Asn147 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2V1/CROC1/UBE2V/UEV1/P/OKcl.19 (His tag) recombinant protein