

CLEC1A antibody



Catalog Number: 109351

Product name

CLEC1A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CLEC1A Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in mouse lung tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of CLEC1A recombinant protein (Accession Number: NM_016511). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109351(CLEC1A antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:109351(CLEC1A antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:109351(CLEC1A antibody) at dilution of 1:100

CLEC1A antibody



Catalog Number: 109351

(under 40x lens)