

FXN antibody



Catalog Number: 110751

Product name

FXN antibody

Specificity

Human; other species not tested.

Antibody description

FXN Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells, human cerebellum tissue, human heart tissue, Jurkat cells, K-562 cells, Raji cells. Observed molecular weight by Western-blot: 14kd

Preparation

This antibody was obtained by immunization of FXN recombinant protein (Accession Number: NM_000144). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:110751(FXN antibody) at dilution of 1:500