

# ACHE Polyclonal Antibody



Catalog Number: 164751

---

**Product name**

ACHE Polyclonal Antibody

**Specificity**

Mouse

**Antibody description**

Polyclonal antibody to ACHE

**Preparation**

Antigen: Recombinant protein of human ACHE

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**


Rabbit IgG

**Applications**

IHC

**Dilutions**

IHC 1:50 - 1:200

**Validations**

Immunohistochemistry - ACHE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse intestine, using ACHE antibody .