

# Anti-GRSF1 Rabbit antibody



Catalog Number: 167751

## Product name

Anti-GRSF1 Rabbit antibody

## Specificity

Human

## Antibody description

GRSF1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human GRSF1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000

IHC: 1/20

ICC/IF: 1/20

FC: 1/50

## Validations



Western blot detection of GRSF1 in HeLa,Ramos cell lysates using GRSF1 Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.