

# CA7 antibody

Catalog Number: 108751

## Product name

CA7 antibody

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CA7 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human stomach tissue. Positive WB detected in human kidney tissue. Observed molecular weight by Western-blot: 30kd

## Validations




## Preparation

This antibody was obtained by immunization of CA7 recombinant protein (Accession Number: NM\_005182). Purification method: Antigen affinity purified.

human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108751(CA7 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108751(CA7 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

# CA7 antibody



Catalog Number: 108751

---

No:108751(CA7 Antibody) at dilution of 1:50 (under 40x lens)