

DCLK1 antibody



Catalog Number: 109751

Product name

DCLK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCLK1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa, 82 kDa

Preparation

This antibody was obtained by immunization of DCLK1 recombinant protein (Accession Number: NM_004734). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109751(DCLK1 Antibody) at dilution of 1:1500