

NDUFA5 antibody



Catalog Number: 113061

Product name

NDUFA5 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

NDUFA5 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, human heart tissue, rat liver tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 13kd




Preparation

This antibody was obtained by immunization of NDUFA5 recombinant protein (Accession Number: NM_005000). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113061(NDUFA5 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:113061(NDUFA5 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human brain using Catalog No:113061(NDUFA5 antibody) at dilution of 1:50 (under 40x lens)