

# Pelgifatamab reference antibody



Catalog Number: 177061

## Product name

Monoclonal

Pelgifatamab reference antibody

## Ig Type

## Specificity

IgG1

Human FOLH1 / PSMA

## Applications

## Antibody description

ELISA

Anti-FOLH1 / PSMA Antibody (Pelgifatamab)

## Validations

## Preparation

Recombinant expression and purified from CHO cells.


## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality



## SDS-PAGE

Anti-FOLH1 / PSMA Antibody (pelgifatamAb) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.