

# Human MGAT2 (His tag) recombinant protein



Catalog Number: 519461

## General Information

### Gene Name Synonym

Alpha-1,6-Mannosyl-Glycoprotein 2-Beta-N-Acetylglucosaminyltransferase, Beta-1,2-N-acetylglucosaminyltransferase II, GlcNAc-T II, NT-II, Mannoside Acetylglucosaminyltransferase 2, N-Glycosyl-Oligosaccharide-Glycoprotein N-Acetylglucosaminyltransferase II, MGAT2

### Protein Construction

A DNA sequence encoding the Human MGAT2 Arg30-Gln447 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 150mM NaCl, pH8.0 .

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human MGAT2 (His tag) recombinant protein