

Integrin alpha-9 antibody



Catalog Number: 111861

Product name

Integrin alpha-9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Integrin alpha-9 Rabbit Polyclonal antibody.
Positive IF detected in A549 cells. Positive WB detected in human liver tissue, A549 cells.
Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of Integrin alpha-9 recombinant protein (Accession Number: BC030198). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111861 (ITGA9 antibody) at dilution of 1:400



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using Catalog No:111861 (Integrin alpha-9 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)