

MRPS9 antibody



Catalog Number: 112861

Product name

MRPS9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MRPS9 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of MRPS9 recombinant protein (Accession Number: NM_182640). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112861(MRPS9 Antibody) at dilution of 1:1000