

# RTN1 (Isoform RTN1-C) antibody



Catalog Number: 114861

## Product name

RTN1 (Isoform RTN1-C) antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RTN1 (Isoform RTN1-C) Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue. Positive IHC detected in human brain tissue, human placenta tissue. Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of RTN1 (Isoform RTN1-C) recombinant protein (Accession Number: NM\_206852). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114861(RTN1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human brain using Catalog No:114861(RTN1 antibody) at dilution of 1:100 (under 40x lens)