

Anti-Beclin 1 Rabbit antibody



Catalog Number: 166861

Product name

Anti-Beclin 1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Beclin 1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Beclin 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Beclin 1 in MCF-7 cell lysates using Beclin 1 Rabbit pAb(1:1000 diluted). Predicted band size:52KDa. Observed band size:52KDa.