

Anti-TPOR/MPL Rabbit antibody



Catalog Number: 168861

Product name

Anti-TPOR/MPL Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TPOR/MPL Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human TPOR

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

1:1000-1:5000

Validations



Western blot detection of TPOR/MPL in K562, C6 cell lysates using TPOR/MPL Rabbit mAb(1:1000 diluted). Predicted band size: 71 kDa. Observed band size: 85 kDa.