

Anti-PRMT7 antibody



Catalog Number: 141161

Product name

Anti-PRMT7 antibody

Specificity

Human PRMT7

Antibody description

Rabbit Polyclonal to Human PRMT7

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human PRMT7, and purified by antigen affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC-P, ICC/IF, IF

Dilutions


IHC-P: 0.1-2 μ g/mL

ICC/IF: 0.5-3 μ g/mL

Validations



Human PRMT7 Immunohistochemistry(IHC)



15704

Human PRMT7 Immunofluorescence(IF) 15705