

CLEC4G antibody



Catalog Number: 107161

Product name

CLEC4G antibody

Specificity

Human; other species not tested.

Antibody description

CLEC4G Mouse Monoclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of CLEC4G recombinant protein (Accession Number: NM_198492). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107161(CLEC4G Antibody) at dilution of 1:1000