

FBXW11 antibody

Catalog Number: 110561

Product name

FBXW11 antibody

IP: 1:200-1:2000

IHC: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FBXW11 Rabbit Polyclonal antibody. Positive IHC detected in human stomach cancer tissue. Positive WB detected in SGC-7901 cells, human stomach tissue, mouse stomach tissue. Positive IP detected in mouse stomach tissue. Observed molecular weight by Western-blot: 61 kDa

Preparation

This antibody was obtained by immunization of FBXW11 recombinant protein (Accession Number: NM_033644). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:110561 (FBXW11 antibody) at dilution of



IP Result of anti-FBXW11 (IP: Catalog No:110561, 4ug; Detection: Catalog No:110561 1:600) with mouse stomach tissue lysate 4000ug.




Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using Catalog No:110561 (FBXW11 antibody).

FBXW11 antibody

Catalog Number: 110561

No:110561(FBXW11 Antibody) at dilution of 1:100
(under 10x lens).



Immunohistochemistry of paraffin-embedded
human stomach cancer tissue slide using Catalog
No:110561(FBXW11 Antibody) at dilution of 1:100
(under 40x lens).