

Cynomolgus OLR1/LOX1 (His Tag) recombinant protein



Catalog Number: 502561

General Information

Protein Construction

A DNA sequence encoding the cynomolgus OLR1 (G7PJT2) (Ser61-Gln273) was expressed with a polyhistidine tag at the N-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant heterodimer of cynomolgus OLR1 comprises 233 amino acids and has a calculated molecular mass of 26.9 KDa. The

apparent molecular mass of the protein is approximately 37 KDa respectively in SDS-PAGE.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus OLR1 / LOX1 Protein (His Tag) SDS-PAGE