

# MCRS1 antibody



Catalog Number: 112561

## Product name

MCRS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MCRS1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, human liver tissue, human spleen tissue, Jurkat cells, mouse spleen tissue, mouse uterus tissue, rat spleen tissue. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of MCRS1 recombinant protein (Accession Number: NM\_006337). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112561(MCRS1 antibody) at dilution of 1:500