

# SPDYA antibody

Catalog Number: 115561

## Product name

SPDYA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SPDYA Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of SPDYA recombinant protein (Accession Number: NM\_001142634). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-SPDYA (IP:Catalog No:115561, 3ug; Detection:Catalog No:115561 1:300) with mouse testis tissue lysate 4000ug.