

# H4N6 HA (His Tag) recombinant protein



Catalog Number: 500961

## General Information

### Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 343) of the influenza A hemagglutinin (A/Swine/Ontario/01911-1/99 (H4N6)) (AAG17429.1), termed as HA1 was expressed, fused with a C-terminal polyhistidine tag.

### Organism

H4N6

### Expression Host

Human Cells

## QC Testing

### Purity

> 97 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Gln 17

### Molecular Mass

The secreted recombinant influenza A H4N6 HA1 subunit (A/Swine/Ontario/01911-1/99 (H4N6)) comprises 338 amino acids and has a predicted

molecular mass of 37.6 kDa. As a result of glycosylation, the apparent molecular mass of the protein is approximately 45-55 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H4N6 (A/Swine/Ontario/01911-1/99) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE