

# Anti-HSP90 Rabbit antibody



Catalog Number: 168961

## Product name

Anti-HSP90 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

HSP90 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Hsp90 alpha

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000

ICC/IF: 1/100

## Validations



Western blot detection of HSP90 in K562,C6,HeLa cell lysates using HSP90 Rabbit mAb(1:1000 diluted).Predicted band size:85kDa.Observed band size:90kDa.