

# Human CASP10/MCH4 (His tag) recombinant protein



Catalog Number: 518961

## General Information

### Gene Name Synonym

Caspase-10, CASP-10, Apoptotic Protease Mch-4, FAS-Associated Death Domain Protein, Interleukin-1B-Converting Enzyme 2, FLICE2, ICE-Like Apoptotic Protease 4, CASP10, MCH4

### Protein Construction

A DNA sequence encoding the Human CASP10/MCH4 Val220-Ile480 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 25mM HEPES, 10mM DTT, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human CASP10/MCH4 (His tag) recombinant protein