

DNASE2B antibody



Catalog Number: 109961

Product name

DNASE2B antibody

Specificity

Human; other species not tested.

Antibody description

DNASE2B Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive IHC detected in human tonsillitis tissue, human prostate hyperplasia tissue. Positive WB detected in Jurkat cells, A549 cells. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of DNASE2B recombinant protein (Accession Number: BC126136). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109961(DNASE2B Antibody) at dilution of 1:2000



Immunofluorescent analysis of A549 cells using Catalog No:109961(DNASE2B Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109961(DNASE2B Antibody) at dilution of 1:50 (under 40x lens)