

# SCN3A Polyclonal Antibody



Catalog Number: 160261

## Product name

SCN3A Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to SCN3A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 274-400 of human SCN3A (NP\_008853.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



## Western blot - SCN3A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SCN3A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.