

# EIF4E3 antibody



Catalog Number: 110261

## Product name

EIF4E3 antibody

## Immunogen

[Human EIF4E3 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EIF4E3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 24 kDa

## Preparation

This antibody was obtained by immunization of EIF4E3 recombinant protein (Accession Number: XM\_017006279). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110261(EIF4E3 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using EIF4E3 antibody Catalog No:110261 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).