

NO66/C14orf169 antibody



Catalog Number: 113261

Product name

NO66/C14orf169 antibody

Specificity

Human; other species not tested.

Antibody description

NO66/C14orf169 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in A431 cells, Jurkat cells. Observed molecular weight by Western-blot: 71 kDa

Preparation

This antibody was obtained by immunization of NO66/C14orf169 recombinant protein (Accession Number: NM_024644). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113261(C14orf169 Antibody) at dilution of



1:600

IP Result of anti-C14orf169 (IP:Catalog No:113261, 4ug; Detection:Catalog No:113261 1:500) with Jurkat cells lysate 4000ug.