

Product name

Anti-Atrx antibody

Dilutions

ICC: 1:200

Specificity

IHC-P: 1:200

Mouse

Validations**Antibody description**

Mouse monoclonal antibody to Atrx

Preparation

This antigen of this antibody was this antibody is produced by immunizing mice with recombinant protein of mouse atrx.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ICC, IHC-P




Fig1: ICC staining ATRX in mouse embryonic cells D3 (red). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2: Immunohistochemical analysis of paraffin-embedded mouse kidney tissue using anti- ATRX mouse monoclonal antibody.