

# Anti-Atrx antibody



Catalog Number: 176661

## Product name

Anti-Atrx antibody

## Specificity

Mouse

## Antibody description

Mouse monoclonal antibody to Atrx

## Preparation

This antigen of this antibody was this antibody is produced by immunizing mice with recombinant protein of mouse atrx.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ICC, IHC-P

## Dilutions

ICC: 1:200

IHC-P: 1:200

## Validations




Fig1: ICC staining ATRX in mouse embryonic cells D3 (red). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2: Immunohistochemical analysis of paraffin-embedded mouse kidney tissue using anti- ATRX mouse monoclonal antibody.